

Views from Lake Placid




A Brief Survey of the Chromatin Landscape

The process of organizing the millions of base pairs of genetic material in the eukaryotic nucleus has profound effects on DNA-dependant events, such as transcription, recombination, replication and repair. DNA is organized by its incorporation into chromatin, the basic subunit of chromatin being the nucleosome. A nucleosome is composed of 147 base pairs of DNA coiled around an octamer of histone proteins, two molecules each of histone H2A, H2B, H3 and H4. Histone H1 associates with chromatin outside the core octamer unit and regulates higher order chromatin structure. Chromatin is packaged into chromosomes by proteins such as the condensins and the chromatin assembly proteins, but chromatin and chromosomes undergo dramatic and dynamic changes in organization in response to a myriad of cellular signals. Chromosomes condense and relax during the cell division process. Damaged DNA adopts a unique structure that facilitates repair. And critical to cell function, most of the genome must remain in a transcriptionally silent state, save for specific combinations of genes which vary significantly between cell types. These processes must be tightly regulated to maintain the integrity of the genome and the proper function of the cell.

Multiple mechanisms exist within the nucleus that allow the function and organization of the genome to be dynamically regulated, responding rapidly to different signals or insults. During the past fifteen years, it has become evident that ATP-dependant enzymes that reorganize chromatin at the level of the nucleosome (Swi/Snf, Rsc, Iswi...) play an important part in the regulation of genomic organization. In mammals, CpG methylation of DNA is also involved in DNA-dependant processes, tied inextricably to other chromatin-based mechanisms of modulating genome function. Most notably though, a large body of evidence dating back forty years or more has accumulated to indicate that post-translational modification of the histone proteins that form the core of the nucleosome is crucial to all genome-based activity (Table 1; see Kouzarides, 2007 for a recent and comprehensive review of the subject). The recruitment of transcriptional regulators and chromosomal proteins brings along enzymes that modify (by either addition or removal) specific functional groups on histones, and these dynamic addition and subtraction events have profound effects on the structure and function of chromatin.

Many of the proteins involved in the deposition or removal of modifications are involved in important cellular process such as cell division, DNA damage repair, cell fate determination, and the maintenance of stem cell pluripotency (see Ringrose and Paro, 2007; Spivakov and Fisher, 2007). So important are these proteins that the alteration in function of a histone modifying enzyme is frequently associated with specific cancers, making histone modifiers attractive targets for therapeutic intervention (see Swamanathan et al., 2007).



Immunodetection of Histone Acetylation in HeLa Cells: Red: Anti-alpha-Tubulin (AM-0191), Green: Anti-acetyl-Lys23 Histone H3 (AR-0104), Blue: DAPI (CK-0156). Far right: Merge.

Generally, acetylation of histones is involved with loosening chromatin structure and the activation of gene expression. Deacetylation is correlated with the tightening of chromatin structure and the repression of expression. The only acetylation mark that appears to have an identified role in specific transcriptional processes is Ac-Lys 16 of histone H4 (see Shia et al, 2006). Histone methylation has a more complex role in controlling access to specific regions of the genome than acetylation. Both lysine and arginine residues in histones can be methylated, and both are subject to multiple forms of methylation (mono-, di-, tri-methylation of lysines; mono-, symmetric or asymmetric di-methylation of arginines). Also, many more examples exist where the methylation of a specific amino acid is associated with a unique chromosomal phenomenon (e.g., lysine 36 methylation of H3 and transcriptional elongation). With regard to histone phosphorylation, many such events are involved in chromosome condensation and segregation. Other marks, such as sumoylation and ubiquitylation have functional roles just now being made clear.

continued...



C. David Allis on [Chromatin Reagents](#)

"Modifications of histones produce a genomic 'bar code' that contains important epigenetic instructions that influence gene expression and cellular development.

The applications employed in studying epigenetics and chromatin biology can be very unforgiving. That's why it's important to use the highest quality reagents."

We agree.

At Lake Placid Biologicals, our chromatin and cell signaling research reagents reach your lab ONLY after rigorous evaluation in ours.

AND, they are backed by unparalleled, real-time technical support and customer service.

Antibodies

Histone Modifications
Chromatin Modifiers
ChIP Validated
IF Optimized

Chromatin IP


ChIP Kits
Protein G-Agarose

Immunostaining

Cell Staining Kits
Custom Services
Fluorescent Secondary
Antibodies

Reagent failure is not an option.

Complete product portfolio inside.



As complex as the histone modification landscape is, it is even more so since the identification of a co-repressor protein complex member, BHC110, as the first histone demethylase. This groundbreaking discovery helped make possible the identification of dozens more demethylating enzymes and has opened up a new field of study for researchers studying chromatin biology and its influence on the regulation of genomic activity (see Shi, 2007). Many of these novel demethylases have been previously identified and studied for their roles in development, or have been found to exist as part of complexes containing widely studied transcription factors. A large amount of work is going regarding demethylases, not the least of which involves identifying the catalytic mechanisms of the different families of demethylases, the substrate specificity and the importance of demethylation in the grand scheme of genome regulation and organization.

As the role of a few of the vast number of histone modifications was uncovered, the existence of a "histone code" was put forth, helping to make sense of the panoply of histone modifications and to try to predict and define cellular events associated with specific modification patterns. The histone code hypothesis suggests that post-translational modifications (or combinations of modifications) confer specific biological functions and that modified histones recruit specialized chromatin-interacting proteins that facilitate the defined function. It is believed that these histone-effector protein interactions give rise to downstream protein recruitment and enzyme and substrate interactions (Table 1). Proteins such as HP1 and Pc bind to repression-associated methylation events on histone H3 to maintain regions of the genome in a repressed state. The Bromo domain binds to acetylated histones and appears to be involved in recruitment of histone modifying complexes. It is likely that other such domains remain to be discovered, which will help identify the molecular mechanisms by which histone modifications exert their influence on the genome.

Table 1:

Histone Variant	Identity vs. Core Histone	Function
H2A.X	96%	DNA damage repair; ser in SQEY motif in mammalian H2AX phosphorylated by ATM and ATR in response to double strand DNA breaks
H2A.Z	60-65%	Transcriptional regulation, DNA damage, chromosome stability
macroH2A	64%	Enriched on inactive mammalian X chromosome; unique N-terminal addition
CENP-A	highly diverged	Centromere-specific; specialized protein for attachment of kinetochore
H3.3	97%	Replaces H3.1 at transcriptionally active regions; varies at 4 amino acids

In addition to the core histones, several unique versions of the histone proteins exist that are not identical to the canonical core proteins. These "variant" histones have unique functions, unique expression patterns and in some cases are restricted to specific sub-nuclear domains (see Table 2 and Henikoff and Ahmad, 2005).

Currently, the chromatin immunoprecipitation (ChIP) technique is the most precise method to identify specific proteins associated with a region of the genome, or conversely, identification of regions of the genome associated with specific proteins. However, much of the recent usage of ChIP has been for identifying and analyzing transcription factor binding sites.


Table 2:

Modification	Associated Pathway	Modification Binding Domain
Acetylation	Transcriptional Regulation; DNA Repair; Chromosome Condensation	Bromo
Methylation-Lysine: mono, di-, tri-	Transcriptional Regulation; DNA Repair; Heterochromatin Formation	Chromo, MBT, Tudor, PHD
Methylation-Arginine: mono-, asym. di-, sym. di-	Transcription	
Phosphorylation: Ser, Thr	Chromosome Condensation; Transcriptional Regulation	14-3-3
Ubiquitylation: Lys	Transcriptional Regulation	
Sumoylation: Lys	Transcriptional Regulation	
ADP-ribosylation: Glu	Transcriptional Regulation	

New techniques that combine the ChIP technique with microarray analysis facilitate the determination of genome-wide distribution patterns for chromosomal proteins and histone modifications. With this combined technique, researchers can ask sweeping questions regarding cellular functions, and place the protein or modification in the larger context of regulatory network interactions. Recent advances in high-throughput DNA sequencing threaten to make ChIP on chip obsolete, as vastly larger and more accurate data sets can be generated (see Barski et al., 2007). As powerful as the ChIP technique is, the chromatin field is in need of other such versatile techniques that can help researchers understand the nature of histone modifications, chromatin organization and their role in biology and disease.

References:

- Barski A et al., (2007) *Cell* **129**:823-827.
 Henikoff S, Ahmad K (2005) *Ann Rev Cell Dev Bio* **21**:133-153.
 Kouzarides T (2007) *Cell* **128**:693-705.
 Ringrose L, Paro R (2007) *Development* **134**:223-232.
 Shi Y, Whetstone JR (2007) *Mol Cell* **25**:1-14.
 Shia WJ et al., (2006) *Genes Dev* **20**:2507-2512.
 Swaminathan V et al., (2007) *Subcell Biochem* **41**:397-428.



Resetting Histone Methylation: Crystal structure of the catalytic domain of the histone demethylase JMJD2A bound to a peptide derived from the amino terminus of histone H3, trimethylated at lysine 9 (yellow carbon atoms). Also shown is the substrate alpha-ketoglutarate (purple) and a Ni(II) cation (red) which is substituted for the Fe(II) normally present. Image courtesy Dr. Ray Trievel, Department of Biological Chemistry, University of Michigan.

Products and Applications

	Description	Applications	Reactivity	Cat. No.	Vial Size	USD
Antibodies and Purified Proteins						
ASH2	Anti-ASH2	WB IF IP	Hu WR*	AR-0175-200	200 µl	\$275
ATM	Anti-p-ATM (Ser1981), Clone 10H11.E12	IF	Hu Ms	AM-0173-200	200 µl	\$275
Copine I	Anti-Copine I	WB	Hu WR*	AR-0167-200	200 µl	\$275
EZH2	Anti-EZH2, polyclonal	WB IF ChIP IHC	Hu Ms	AR-0163-200	200 µl	\$275
Hbo1	NEW Anti-Hbo1, polyclonal	WB	Hu Ve*	AR-0197-200	200 µl	\$275
HDAC	Anti-HDAC1, clone 10E2	WB IF IHC IP ChIP	Hu Ms Rt	AM-0131-200	200 µl	\$275
	Anti-HDAC2, clone 3F3	WB IF IHC IP ChIP	Hu Ms Rt	AM-0132-200	200 µl	\$275
Histone H2A	Anti-Histone H2A, acidic patch	WB	Hu WR*	AR-0189-200	200 µl	\$275
	NEW Anti-phospho-Ser1 Histone H2A / Histone H4	WB	Hu WR*	AR-0200-200	200 µl	\$275
	NEW Anti-acetyl-Lys5 Histone H2A	WB	Hu WR*	AR-0201-200	200 µl	\$275
	NEW Anti-acetyl-Lys9 Histone H2A	WB	Hu WR*	AR-0256-200	200 µl	\$275
Histone H2AX	Anti-phospho-Ser139 Histone H2AX	WB IF	Hu Ma*	AR-0149-200	200 µl	\$275
Histone H2A.Z	NEW Anti-Histone H2A.Z	WB IF	Hu Ma*	AR-0255-200	200 µl	\$275
Histone H2B	Anti-Histone H2B	WB IF	Hu Ma*	AR-0139-200	200 µl	\$275
	NEW Anti-acetyl-Lys5 Histone H2B	WB IF	Hu WR*	AR-0228-200	200 µl	\$275
	NEW Anti-acetyl-Lys16 Histone H2B	WB IF	Hu WR*	AR-0227-200	200 µl	\$275
Histone H3	Anti-Histone H3, CT	WB IF ChIP	Hu Sc WR*	AR-0144-200	200 µl	\$275
Histone H3 acetyl	Anti-acetyl Histone H3	WB	Hu WR*	AR-0143-200	200 µl	\$275
	Anti-acetyl-Lys9 Histone H3	WB	Hu WR*	AR-0102-200	200 µl	\$275
	Anti-acetyl-Lys14 Histone H3	WB	Hu WR*	AR-0142-200	200 µl	\$275
	Anti-acetyl-Lys18 Histone H3	WB IF	Hu WR*	AR-0103-200	200 µl	\$275
	Anti-acetyl-Lys23 Histone H3	WB IF	Hu WR*	AR-0104-200	200 µl	\$275
	Anti-acetyl-Lys27 Histone H3	WB ChIP	Hu Sc WR*	AR-0105-200	200 µl	\$275
	Anti-acetyl-Lys27 Histone H3	WB IF ChIP	Hu Sc WR*	AR-0129-50	50 µl	\$150
Histone H3 phospho	Anti-phospho-Ser10,28 Histone H3	WB	Hu WR*	AR-0140-200	200 µl	\$275
	NEW Anti-phospho-Thr11 Histone H3	WB	Hu WR*	AR-0257-200	200 µl	\$275
	NEW Anti-phospho-Ser28 Histone H3	WB	Hu WR*	AR-0198-200	200 µl	\$275
Histone H3 methyl	NEW Anti-dimethyl-Lys4 Histone H3	WB IF	Hu WR*	AR-0168-200	200 µl	\$275
	NEW Anti-trimethyl-Lys4 Histone H3	WB IF ChIP	Hu WR*	AR-0169-200	200 µl	\$275
	Anti-dimethyl-Lys9 Histone H3	WB IF ChIP	Hu WR*	AR-0108-50	50 µl	\$150
	NEW Anti-trimethyl-Lys9 Histone H3	WB IF	Hu WR*	AR-0170-200	200 µl	\$275
	Anti-trimethyl-Lys27 Histone H3	WB IF	Hu WR*	AM-0150-200	200 µl	\$275
	NEW Anti-trimethyl-Lys27 Histone H3	WB ChIP	Hu WR*	AM-0174-200	200 µl	\$275
	NEW Anti-trimethyl-Lys27 Histone H3	WB IF	Hu WR*	AR-0171-200	200 µl	\$275
	NEW Anti-trimethyl-Lys27 Histone H3	WB IF	Hu WR*	AR-0199-200	200 µl	\$275
	NEW Anti-monomethyl-Lys79 Histone H3	WB IF ChIP	Hu WR*	AR-0172-200	200 µl	\$275
	NEW Anti-dimethyl-Lys79 Histone H3	WB	Hu WR*	AR-0177-200	200 µl	\$275
Histone H4 acetyl	Anti-tetra-acetyl Histone H4	WB ChIP	Hu Sc WR*	AR-0109-50	50 µl	\$150
	Anti-tetra-acetyl Histone H4	WB IF	Hu Sc WR*	AR-0120-50	50 µl	\$150
	Anti-acetyl-Lys5 Histone H4	WB IF ChIP	Hu Sc WR*	AR-0119-200	200 µl	\$275
	Anti-acetyl-Lys8 Histone H4	WB	Hu WR*	AR-0138-200	200 µl	\$275
	Anti-acetyl-Lys12 Histone H4	WB IF ChIP	Hu Sc WR*	AR-0106-200	200 µl	\$275
	Anti-acetyl-Lys16 Histone H4	WB ChIP	Hu Sc WR*	AR-0107-200	200 µl	\$275
Histone H4 methyl	Anti-monomethyl-Lys20 Histone H4	WB IF	Hu WR*	AR-0134-200	200 µl	\$275
	Anti-dimethyl-Lys20 Histone H4	WB IF	Hu WR*	AR-0135-200	200 µl	\$275
	Anti-dimethyl-Lys20 Histone H4	WB	Hu WR*	AM-0165-200	200 µl	\$275
	Anti-trimethyl-Lys20 Histone H4	WB IF ChIP	Hu WR*	AR-0136-200	200 µl	\$275
L3MBT	Anti-L3MBTL1	WB IP ChIP	Hu Ve*	AR-0160-200	200 µl	\$275
LexA	Anti-LexA DNA Binding Domain	WB IP	tagged fusions	AR-0166-200	200 µl	\$275
LSD1	Anti-LSD1	WB IP	Hu Ve*	AR-0162-200	200 µl	\$275

continued...

Antibodies (AR-#### are rabbit polyclonal; AM-#### are mouse monoclonal)

Applications Key: **WB**-Western **IP**-Immunoprecipitation **IHC**-Immunohistochemistry
IF-Immunofluorescence **ChIP**-Chromatin IP **FC**-Flow cytometrySpecies Key: **Hu**-Human **Ms**-Mouse **Rt**-Rat **Sc**-Budding yeast **Sp**-Fission yeast
Dr-Drosophila **Ve**-Vertebrate **Ma**-Mammal **Eu**-Eukaryote *-predicted**Reagent failure is not an option.**

Additional vial sizes are available for most antibodies (including 5 µl evaluation samples for \$25).

Products and Applications

	Description	Applications	Reactivity	Cat. No.	Vial Size	USD
MBP	Myelin Basic Protein, dephosphorylated	Kinase Assay		PN-0151-10	10 mg	\$275
MeCP2	NEW Anti-MeCP2	WB	Hu Ms	AR-0195-10	200 µl	\$275
MEP50	Anti-MEP50	WB	Hu Ma*	AR-0145-200	200 µl	\$275
Nucleolin	NEW Anti-Nucleolin	WB IF	Hu	AM-0196-200	200 µl	\$275
PP2A	NEW Anti-PP2A	WB	Hu WR*	AR-0190-200	200 µl	\$275
Rad51	Anti-Rad51	WB IF IHC IP	Hu	AR-0100-200	200 µl	\$275
Rad52	Anti-Rad52	IF IHC	Hu	AR-0101-200	200 µl	\$275
RbAp46/48	NEW Anti-RbAp46/48	WB ChIP	Hu WR*	AR-0178-200	200 µl	\$275
SNF2	Anti-SNF2h	WB IP	Hu Ma	AM-0152-200	200 µl	\$275
Tubulin	NEW Anti-alpha-Tubulin, clone 5-B-1-2	WB IF IP	Hu Ms WR*	AM-0191-200	200 µl	\$275
Chromatin Immunoprecipitation (ChIP) Products						
ChIP	ChIP Protein G-Agarose			MX-0110	2.5 ml	\$100
	Chromatin Immunoprecipitation Kit			KA-0137-25	25 assays	\$185
	Chromatin Immunoprecipitation Kit			KA-013-100	100 assays	\$500
MAX Stain™ Products for Immunofluorescence						
MAXfluor™ Mounting Media	NEW MAXfluor™ for Fluorescence	three formulations designed for optimizing immunofluorescence staining results		SA-0148	2 ml	\$75
	NEW MAXfluor™ with DAPI			SA-0156	2 ml	\$75
	NEW MAXfluor™ with Propium Iodide			SA-0157	2 ml	\$75
	NEW MAXfluor™ Trio Kit (SA-148, -156, -157)			KP-0158	Kit	\$200
MAXblock™	NEW MAXblock™ Blocking Medium	most effective blocking reagent available		SA-0146	150 ml	\$85
MAXbind™	NEW MAXbind™ Staining Medium	staining reagent for IF optimization		SA-0147	250 ml	\$150
MAXwash™	NEW MAXwash™ Immunostain Washing Medium	high efficiency washing		SA-0226	1 L	\$75
MAXpack™	NEW MAXpack™ Immunostain Medium Kit	includes MAXblock™, MAXbind™ & MAXwash™		SA-0148	Kit	\$285
MAX Stain™ Slip and Slide	NEW MAX Stain™ Slip and Slide	50 ea. slides and #1.5 coverslips		KP-0188	50 ea.	\$80
	NEW MAX Stain™ Slip and Slide, acid etched	50 ea. slides and #1.5 coverslips		KP-0221	50 ea.	\$100
	NEW MAX Stain™ Slip and Slide, Poly-D-Lysine	50 ea. slides and #1.5 coverslips		KP-0222	50 ea.	\$140
MAX Stain™ Deep Space Probes	NEW Orion Green 488 Anti-Mouse IgG	maximally cross-adsorbed goat anti-mouse		CA-0205	250 µl	\$120
	NEW Orion Green 488 Anti-Rabbit IgG	maximally cross-adsorbed goat anti-rabbit		CA-0206	250 µl	\$120
	NEW Triskelion Yellow 549 Anti-Mouse IgG	maximally cross-adsorbed goat anti-mouse		CA-0208	250 µl	\$120
	NEW Triskelion Yellow 549 Anti-Rabbit IgG	maximally cross-adsorbed goat anti-rabbit		CA-0209	250 µl	\$120
	NEW Regulus Red 594 Anti-Mouse IgG	maximally cross-adsorbed goat anti-mouse		CA-0251	250 µl	\$120
	NEW Regulus Red 594 Anti-Rabbit IgG	maximally cross-adsorbed goat anti-rabbit		CA-0252	250 µl	\$120
	NEW Vulcan Red 649 Anti-Mouse IgG	maximally cross-adsorbed goat anti-mouse		CA-0230	250 µl	\$120
	NEW Vulcan Red 649 Anti-Rabbit IgG	maximally cross-adsorbed goat anti-rabbit		CA-0231	250 µl	\$120
	NEW Feris VI Far Red 680 Anti-Mouse IgG	maximally cross-adsorbed goat anti-mouse		CA-0233	250 µl	\$120
	NEW Feris VI Far Red 680 Anti-Rabbit IgG	maximally cross-adsorbed goat anti-rabbit		CA-0234	250 µl	\$120
	NEW Ultra Far Red 800 Anti-Mouse IgG	maximally cross-adsorbed goat anti-mouse		CA-0236	250 µl	\$120
	NEW Ultra Far Red 800 Anti-Rabbit IgG	maximally cross-adsorbed goat anti-rabbit		CA-0237	250 µl	\$120
MAX Stain™ CytoMark	NEW Orion Green Anti-alpha Tubulin	IF direct conjugate antibody - mouse IgG		CA-0214	100 µl	\$250
	NEW Orion Green Anti-Trimethyl Lys27 Histone H3	IF direct conjugate antibody - rabbit IgG		CA-0215	100 µl	\$250
MAX Stain™ IF Kit	NEW MAX Stain™ Universal Immunofluorescence Kit	includes 2 6-well plates with coverslips 2 Deep Space Probes (your choice), MAXblock™ MAXbind™ staining medium, MAXwash™ washing medium and MAXfluor™ mounting medium		KA-0179	1 kit	\$255

Ordering Information

Web Store: www.lpbio.com (accepts credit cards and purchase orders)
Phone: 1-866-572-4670 or 518-523-7025
Fax: 518-523-7029
email: orders@lpbio.com



1915 Saranac Ave
 Lake Placid, NY 12946 USA
 www.lpbio.com

www.lpbio.com

1-866-572-4670

518-523-7025